

ceraMotion® Zr



Photo: © Christian Ferrario

ceraMotion® Zr Paste Glaze

ceraMotion® Zr glaze material in paste form



D
DENTAURUM

ceraMotion® Zr Paste Glaze

Glaze material for zirconia and lithium disilicate frameworks

Paste Glaze

Shade code: PGL transpa, 3 g, REF 250-271-03

Advantages

- ▷ ready to use
- ▷ fluorescent
- ▷ reliable wetting of the material surface
- ▷ targeted and controlled application possible
- ▷ only one glaze firing required
(homogeneous degree of glaze already after the first glaze firing)
- ▷ particularly suitable for monolithic works made of zirconia and lithium disilicate
- ▷ uniform degree of glaze even when staining and build-up techniques are combined
- ▷ easy handling

Firing table

	Start temp. (°C/°F)	Drying time (min)	Heat rate (°C/min)	Vacuum start (°C/°F)	Vacuum end (°C/°F)	Firing temp. (°C/°F)	Holding time (min)
Glaze firing*	500 / 932	6	55	500 / 932	750 / 1382	750 / 1382	1 without vacuum

* glaze firing can be completed with or without vacuum
Please observe the instructions provided by the manufacturer when fabricating the framework.



„To me, ceraMotion® means purity and art in designing dental aesthetics.“

Master Dental Technician Germano Rossi, Italy



Photo: ©: Germano Rossi



Minimal effort –
enhanced aesthetics



Dentaurum Group

Germany | Benelux | España | France | Italia | Switzerland | Australia | Canada | USA
and in more than 130 countries worldwide.



DENTAURUM
QUALITY
WORLDWIDE
UNIQUE

➔ For more information on our products and services, please visit www.dentaurum.com

Date of information: 04/16

Subject to modifications



www.dentaurum.com



Like us on Facebook!



Visit us on YouTube!

D
DENTAURUM