Multi 800

manual


Model : MULTI-800 / MULTI-800S

※ Read this User's Manual thoroughly for safety before preparing the for use.


1. Precaution

- Install the unit at the working place and do not use it in a dirty or very hot and humid place. (-20℃ ~ 40℃, 0 ~ 90% RH)
- Please be careful not to drop the handpiece. It will damage ball bearing or weaken durability of motor.
- Please do not allow water or other liquids to spill onto and into the handpiece and control box.
- After replacing a bur, please make sure to close the chuck handle. Also, do not to open the chuck handle, while the handpiece is operating.
- While the handpiece is not in operation, it is recommended to keep the bur into the collet chuck after cleaning up.

2. Component parts


3. Part Introduction

4. Main Power

- 1) Check the electric rating label(230V~, 50Hz) affixed to the rear of the unit.
- 2) Use only the fuse with allowable electric voltage (1.6AL 250V) for the Fuse Holder


- Never plug or unplug the Power Cord with wet hands to avoid electric shock.
- Be sure to use only the power earthed.

5. Installation

5-1. Control unit


When using the standard type as a desktop, install it on the work table.
When using it as a knee control unit, install it on the lower wall of the work table at knee height.

- The Control Unit should rest on a flat, steady surface
- Use the Control Unit in locations that have indoor temperature(-20℃ ~ 40℃)
Avoid using the Control Unit in an excessively dusty, warm or humid location, due to electric shock or damage to the unit.

5-2. Connection

- 1) Insert the pedal cord plug firmly after confirming the direction of the foot pedal connector on the back of the control unit.
- 2) Insert the motor cord plug firmly after confirming the direction of the motor connector on the back of the control unit.
- 3) Insert the power cord plug firmly into the power connector on the back of the unit.

- Make sure that all cords are connected safely and properly before plugging

5-3. Bur(Bit) insertion and removal

- 1) To remove a bur, please turn a chuck handle clockwise until it clicks to stop, and then pull out the bur.
- 2) To insert the bur, please turn the chuck handle counterclockwise until it clicks to stop.
- 3) Please be sure not to use bent, asymmetrical or damaged burs. Always use standard burs.

- For 40,000 RPM or more, it is required to use the recommended bur only.
- Always remove dirt before replacing a test bur/polishing point.
- After replacing a test bur/polishing point, it is necessary to make sure that the Collet Chuck is securely to make sure that the Collet Chuck is securely held before running the motor.
- While the motor is running, DO NOT try to adjust the Chuck Handle to prevent the unit from being damaged.
- When the Motor Handpiece is not in use, it is recommended for a test bur/polishing point to be kept inserted into the Motor Handpiece.
- Always put the Motor Handpiece on its stand while in use.
Special care must be taken not to drop it on the floor

6. Operation

6-1. Button setting

Rotation direction setting

- Press button to select the rotation direction.
- Forward : GREEN LED / Reverse : ORANGE LED

Hand / Foot mode setting

- Press the button to switch to hand/foot (knee) mode.
- Hand mode : LED OFF / Foot (knee) mode: GREEN LED

Handpiece operation /Auto cruise

- Hand mode : Press the button to operate or stop the handpiece.
Handpiece operation: GREEN LED / Handpiece stop: LED OFF
- Foot mode : When foot (knee) pedal is maintained at a constant speed for 2 seconds after setting foot mode, it operates in auto cruise mode. Auto cruise mode is disengaged by pressing the pedal again or pressing the button again.
Auto cruise mode setting: ORANGE LED / Auto cruise mode release:GREEN LED

6-2. Operating procedure

- 1) Connect the power cord plug to a power outlet.
- 2) Set the speed control dial to the lowest speed position.
- 3) Turn on the power switch.
- 4) If you turn the speed control dial on the front of the control box, the speed display window turns on.
- 5) Set the mode with the hand/foot mode button.
- 6) If you press the foot (knee) pedal in foot mode, the motor rotates.
- 7) In hand mode, set the speed control dial to the desired speed and press the handpiece operation button to rotate the motor.

- Wear safety glasses and mask whenever doing cutting or grinding works.
- Do not exceed the maximum operating speed of the Control Unit because its operation at too high speed may result in damage to the unit, as well as personal injury.
- Special care must be taken not to drop the Motor Handpiece or the Control Unit on the floor.
- Before turning off the Power Switch, it is surely required to make sure that the motor has been stopped.
- Unplug the Power Cord when the unit is not in use.


7. Safety protection system

- When the motor stops due to overload, the motor rotates only when the speed control dial is turned to the lowest speed and then turned to the desired speed again.

8. Maintenance


8-1. How to mount and dismount the Collet Chuck

- Turn the Collet Chuck counterclockwise while opened to dismount. If the Collet Chuck can not be turned manually use the Chuck Wrench to dismount it.
- To mount the Collet Chuck, it is required to first open the Chuck Handle and then insert the Collet Chuck and the Bur

8-2. How to disassemble the Chuck Handle and the Motor

- The Chuck Handle and the Motor are disassemble by continuously turning the Chuck Handle counter clockwise. (As much as possible, please ask an expert for separation.)

8-3. How to disassemble the Motor and the Motor Cord

- The Motor Cord is disassembled by continuously turning the Motor counterclockwise.

9. Error code

Error code	Status	Cause
"E1"	Motor sensor error	<ul style="list-style-type: none"> • The Collet Chuck is opened. • The Motor Handpiece is mechanically defective.
"E2"	Overcurrent detection error	• handpiece overload
"E3"	Overvoltage Detection error	• The circuit of the Control Unit is defective.
"E4"	Overheat error	• Overheating inside the control unit
"E5"	Power supply voltage abnormal	• Power supply voltage detection error at the output terminal of the transformer
"E6"	Handpiece connection error	• The motor cable is missing or has poor contact
"E7"	Abnormal operation error	• Overload or sudden stop after normal operation of the motor

10. Troubleshooting

Troubles	Check points	Corrective actions
The Pilot Lamp Does not light.	Make sure that the Power Cord is correctly plugged.	Correctly plug the Power Cord.
	Check that the fuse is blown.	Replace with a recommended fuse. If the fuse is still blown again, then repair it.
	Make sure that the Power Switch is properly working.	Repair if the Power Switch is troubled.
The Motor does not run.	"E1" Check that the Collet Chuck is opened.	If opened, close the Collet Chuck. If the same error is still displayed, then repair it.
	"E1" Check if the shaft of the Motor Handpiece is rotated by fingers.	If not rotated while closed, the spindle of the Motor Handpiece is troubled, Then, repair it.
	"E2" Check if the Handpiece has been excessively used under the condition of too high load.	Press the Auto cruise button or restart running The Control Unit after turning off the Power Switch. If the same error is still displayed, then repair it..
	"E3" The circuit of the Control Unit it is defective.	Repair it.
	"E4" Check if indoor temperature is too warm in a work place with the control unit.	Try resting for 10 minutes in a cool place and then restarting. If the same error code continues to appear, request repair.
	"E5" Power supply voltage abnormal	Check the power supply voltage of the controller and use it according to the input voltage.
	"E6" Check that the Motor Cord is plugged in proper condition.	Correctly plug the Motor Cord. If the same error is Continuously displayed, replace the Motor cord.
"E7" Abnormal operation error	After turning off the controller, wait 5 minutes and try again.	

The Foot Pedal does not work.	Make sure that the Foot Pedal Cord is properly plugged	Plug the Foot Pedal Cord correctly.
	The Hand/Foot Selector Switch is set by "HAND"	Set the Hand/Foot Selector Switch to "FOOT"
The Motor Handpiece is overheated while in operation.	The ball bearings are defective.	Repair the Motor Handpiece.
The Motor handpiece is excessively vibrated and with serious noise.	Make sure that the Collet Chuck is free from dirt.	Thoroughly clean the Collet Chuck after disassembled.
	Check if the specified test bur/polishing Point is used.	For high speed, use the specified test bur/Polishing point only.
	Check if the test bur/polishing Point is bent.	Replace with new one.
The test bur/polishing point is released.	Make sure that the ball bearing is in good condition.	Repair the Motor Handpiece.
	Check if the Collet Chuck is loose.	If loose, tighten the Collet Chuck.

11. Specifications

11-1 Control Unit

Model	Input	Weight (g)	Dimension(mm)		
			W	D	H
MULTI-800 (Knee control type)	100V - 120V~, 50/60Hz 220V - 240V~, 50/60Hz	2,550	112.5	237	218
MULTI-800S (Standard type)	100V - 120V~, 50/60Hz 220V - 240V~, 50/60Hz	2,230	82.5	225	202.5

11-2 Handpiece


Model	Weight (g)	Dimension (mm)	
		L	Ø
SDE-BS50C	215	160	27

11-3 Foot Pedal

Model	Weight (g)	Dimension(mm)		
		W	D	H
SDE-FS80N	590	117	171	62

12. Allowable rotation speed

- This chart is a general idea of the maximum rotation speed of conventional burs. The allowable rotation speed can be slight different depending on bur materials or types.

※ Notice

- 1) Make sure not to exceed the maximum rotation speed specified by the manufacturer or distributor for the bur being used.
- 2) When the rotation speed of a large bur, even though the final speed is within the allowable range, it may bend or damage the bur and it can cause an accident.
- 3) Never use bent, asymmetrical or damaged burs.
- 4) Always use standard burs.
- 5) Be sure to insert the bur shank **all** the way, until it touches the back of the chuck.
- 6) There are variations in bur shank diameter, depending on the particular bur used. Before using a bur, be sure to adjust the chuck.

13. Warranty

- 1 year
- We do not guarantee its quality in case of users' carelessness.

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