## R-020202

E.....	MUFLA DE VERTIDO PARA SILICONA O GELATINA.....	4
GB.....	POURING FLASK FOR SILICONE/GEL.....	9
F.....	MOUFLE DE DUPLICATION POUR GELATINE.....	14
D.....	GUSSMUFFEL FÜR SILIKON ODER GELATINE.....	19
P.....	MUFLA PARA SILICONE OU GEL.....	24
I.....	MUFFOLA DA COLATA PER SILICONE O GELATINA.....	29

**MESTRA®**

TALLERES MESTRAITUA S.L.

Txori-erri Etorbidea, 60

Tel. (+34)944530388 - Fax (+34)944711725

E-mail: [mestra@mestra.es](mailto:mestra@mestra.es) - [www.mestra.es](http://www.mestra.es)

48150 SONDIKA - BILBAO - ESPAÑA

Recambios disponibles / **Spare parts available** / *Pièces détachées* / Lieferbare Ersatzteile / **Peças sobressalentes disponíveis** / *Ricambi disponibili*

020202-05	Pieza con tres tapones (1 U.) / <b>Three-plug piece (1 u.)</b> / <i>Pièce avec trois bouchons (1 U.)</i> / Teil mit drei Stopfen (1 St.) / <b>Peça com três tampas (1 U.)</b> / <i>Pezzo con tre tappi (1 U.)</i>
050150-02	Imán (1 U.) / <b>Magnet (1 u.)</b> / <i>Aimant (1 U.)</i> / Magnet (1 St.) / <b>Íman (1 U.)</b> / <i>Magnete (1 U.)</i>
050150-04	Arandelas para imán (10 U.) / <b>Magnet washers (10 u.)</b> / <i>Rondelles pour aimant (10 U.)</i> / Haftscheiben für Magneten (10 St.) / <b>Arandelas para íman (10 U.)</b> / <i>Rondelle per magnete (10 U.)</i>
020202-04	Gomas cortas (bolsa 50 U.) / <b>Short rubber bands (50 u. bag)</b> / <i>Elastiques courts (bolsa 50 U.)</i> / Kurze Gummis (Tüte 50 St.) / <b>Borrachas pequenas (pacote 50 U.)</b> / <i>Elastici corti (sacchetto 50 U.)</i>
020202-08	Gomas largas (bolsa 50 U.) / <b>Long rubber bands (50 u. bag)</b> / <i>Elastiques longs (bolsa 50 U.)</i> / Lange Gummis (Tüte 50 St.) / <b>Borrachas grandes (pacote 50 U.)</b> / <i>Elastici lunghi (sacchetto 50 U.)</i>
020202-06	Tubo de corte para bebederos Ø 5 mm (1 U.) / <b>Tube for hole cutting Ø 5 mm (1 u.)</b> / <i>Tubes pour la confection des canaux Ø 5 mm (1 U.)</i> / Schneidrohr für Eingüsse Ø 5 mm (1 St.) / <b>Tubo de corte para canais entrada do molde Ø 5 mm (1 U.)</b> / <i>Tubo di taglio per entrate Ø 5 mm (1 U.)</i>
020202-07	Tubo de corte para bebederos Ø 7 mm (1 U.) / <b>Tube for hole cutting Ø 7 mm (1 u.)</b> / <i>Tube pour la confection des canaux Ø 7 mm (1 U.)</i> / Schneidrohr für Eingüsse Ø 7 mm (1 St.) / <b>Tubo de corte para canais entrada do molde Ø 7 mm (1 U.)</b> / <i>Tubo di taglio per entrate Ø 7 mm (1 U.)</i>

- A Pieza con tres tapones / **Three-plug piece** / *Pièce avec trois bouchons* / Peça com três tampas / **Pezzo con tre tappi** / *Teil mit drei Stopfen*
- B Tapa superior de vertido / **Top pouring cover** / *Partie supérieure du moufle* / Tampa superior de despejo / **Coperchio superiore di colata** / *Oberer Gießdeckel*
- C Enganche para goma de cierre de la mufla / **Hook for rubber band to close the flash** / *Crochet d'arrimage des élastiques de fermeture du moufle* / Engate para borracha de fecho da mufla / **Attacco per elastico di chiusura della muffola** / *Haken für Schließgummi der Muffel*
- D Imán / **Magnet** / *Aimant* / Íman / **Magnete** / *Magnet*
- E Enganche para goma de cierre de caja / **Hook for rubber band to close the box** / *Crochet d'arrimage de l'élastique pour la fermeture de la boîte.* / Engate para borracha de fecho da caixa / **Attacco per elastico di chiusura scatola** / *Haken für Schließgummi der Dose*
- F Caja de limpieza / **Cleaning box** / *Boîte de nettoyage des dents et base du moufle* / Caixa de limpeza / **Scatola di pulizia** / *Reinigungsdose*
- G Pestaña para apertura de la mufla / **Tab for opening the flash** / *Bec pour l'ouverture du moufle* / Cavilha para a abertura da mufla / **Flangia per apertura della muffola** / *Flansch zum Öffnen der Muffel*
- H Pestaña para apertura caja de limpieza / **Tab for opening the cleaning box** / *Bec pour l'ouverture de la boîte de nettoyage* / Cavilha para a abertura da caixa de limpeza / **Flangia per apertura scatola di pulizia** / *Flansch zum Öffnen der Reinigungsdose*
- I Arandelas para imán (5 U.) / **Magnet washers (5 u.)** / *Rondelles pour aimant (5 U.)* / Arandelas para íman (5 U.) / **Rondelle per magnete (5 U.)** / *Haftscheiben für Magneten (5 St.)*
- J Gomas cortas (2 U.) / **Short rubber bands (2 u.)** / *Elastiques courts (2 U.)* / Borrachas pequenas (2 U.) / **Elastici corti (2 U.)** / *Kurze Gummis (2 St.)*
- K Gomas largas (2 U.) / **Long rubber bands (2 u.)** / *Elastiques longs (2 U.)* / Borrachas grandes (2 U.) / **Elastici lunghi (2 U.)** / *Lange Gummis (2 St.)*
- L Pestaña para apertura de los tapones / **Tab for opening the plugs** / *Bec pour le deboîtage des bouchons* / Cavilha para abertura das tampas / **Flangia per apertura dei tappi** / *Flansch zum Öffnen der Stopfen*
- M Tubos de corte para bebederos / **Tubes to cut the holes for resin pouring** / *Tubes pour la confection des canaux de coulée de la résine* / Tubos de corte para canais entrada do molde / **Tubi di taglio per entrate** / *Schneidrohr für Eingüsse*

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Fig. 14

The R-020202 duplication flask for silicone or gel has been specially designed for the “resin pouring” technique. A single device that features all the elements required to perform this well-known technique step by step:

- Worktable with central magnet.
- Flask for silicone/gel pouring and filling.
- Side holes for resin pouring.
- Plugs for cutting the holes for resin pouring.
- Pierced box for steam cleaning the dental pieces.
- Tubes to cut the holes for resin pouring.

The flask is supplied with a magnet, magnet washers, rubber bands for closing and two stainless steel tubes to cut the holes for resin pouring.

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### 2. Preparing the gel/silicone mould:

- 2.1. Place the model on the magnet in the base (D) (Fig. 1).
- 2.2. Put the top cover (B) into place and fit the long rubber band to the hooks (C) (Fig. 2), to make sure it will not open. Place the three-plug piece (A) with the arrow upwards.

Fig. 1


Fig. 2

2.3. Prepare the silicon or gel mixture.

2.4. Pour the gel or silicone through the top central hole until it starts to overflow through the other three holes (Fig. 3). This way a better retention will be achieved.

2.5. Leave the flask to rest for the product to harden.

Fig. 3

2.6. Once the product has hardened, release the rubber band, remove the lower part of the flask by pulling the tab (G) and remove also the three-plug piece by pulling the tab (L).

2.7. Remove the model and leave the silicone/gel mould inside the flask (Fig. 4).

Fig. 4

### 3. Preparing the dental pieces:

- 3.1. Open the cleaning box (F) by pulling the tab (H). Remove the dental pieces from the model one by one and place them in order in each of the spaces in the base (Fig. 5).


Fig. 5


Fig. 6

- 3.2. Close the cleaning box. Fit the short rubber band into the hooks (E) to prevent the box opening and deforming the lid. Try to cover the holes as little as possible (Fig. 6).
- 3.3. Clean the dental pieces in the washing machine or steam clean them through the multiple holes. Clean also the wax remaining on the model (Fig. 7).


Fig. 7

- 3.4. When the dental pieces are clean and dry, insert each of them in its place in the silicone/gel inside the flask (Fig. 8).

Fig. 8

Fig. 9

- 3.5. Use the tubes (M) to cut the three holes for pouring the resin (Fig. 9).

#### 4. Preparing the model:

- 4.1. Make sure the plaster model is completely clean and all the wax has been removed.
- 4.2. Apply a mould-releaser to the plaster model with a paintbrush (Fig. 10)

Fig. 10

Fig. 11

- 4.3. Put the plaster model into its place in the flask (Fig. 11).

## 5. Obtaining the resin piece

- 5.1. Close the flask and fit the long rubber band into the hooks (C) to make sure it will not open.
- 5.2. Prepare the resin and pour it through one of the pouring holes until the resin overflows through the other two holes (Fig. 12).

Fig. 12

Fig. 13

- 5.3. Put the flask inside the pressure polymerizer with the holes for pouring the resin upwards. Make sure the water does not reach the resin (Fig.13). Follow the directions about pressure and temperature given by the resin manufacturer.
- 5.4. When the polymerization is over, open the flask, empty it and cut the silicone to remove the piece (Fig. 14).

Fig. 14

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