Mode d'emploi

Sebond pink

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Sebond pink

Product description:

Sebond pink is a light-curing paste opaquer for covering model casts with an integral adhesive component.

Use of an additional bonding system is therefore unnecessary.

Use:

The framework is worked as usual after casting. The surface is sand-blasted with aluminium oxide 50-110 µm, 2 bar, and blown off with oil-free pressurised air.

Sebond pink is applied thinly with a firm short-haired brush. At least two coats are needed to obtain complete coverage. A non-covering opaquer layer has a negative effect on the shade. The flow characteristics of the opaquer can be optimised by mixing it on a mixing block. Ensure that it is applied particularly thinly in the region of the retentions. Each layer polymerises separately.

Spektra 2000: 3 min.
Spektra LED: 1 min.

Sebond pink is available in the shades light and dark. The two colours can be mixed together to obtain the individual shade.

Sources of error / avoiding errors:

The surfaces of the framework must be cleansed of grease layers and polish residues. The opaquer must be applied only in very thin layers. Sebond pink cannot be mixed with other opaquer liquids or powders. Do not use any intermediate polymerising devices for polymerisation.

Opaque layers applied too thickly prevent optimum polymerisation and thus weaken the bonding. Prolonging the polymerisation time does not produce a greater depth of curing.

Working time of the mass:

1-3 min depending on the light.

Storage:

The material should be stored at 10-25 °C. Close the syringe well immediately after use and protect from direct light.

Shelf life:

The maximum shelf life is printed on the label of each syringe. Do not use after the expiry date. Syringe contents: 5 g

Guarantee:

Our technical instructions for use, whether provided orally, in writing or by practical instruction, are based on our own experience and can therefore be regarded only as guidelines. Our products undergo continuous development. We therefore reserve the right to make alterations.

Side effects:

Undesirable side effects of this medicinal product can be expected extremely rarely if it is handled and used correctly. However, immune reactions (e.g. allergy) or local discomfort cannot be completely excluded. If you become aware of undesirable side effects, even in doubtful cases, please notify us.

Contraindications / interactions:

This product must not be used or must be used only under the strict supervision of the treating doctor or dentist in the event of patient hypersensitivity to one of the constituents. In these cases, the composition of the medicinal product supplied by us is available on request. Known cross reactions or interactions between the medicinal product and other materials already present in the mouth must be taken into account by the dentist when using it.

Note:

Please give all the above information to the treating dentist if you are employing this medicinal product for a special model. When using it, please follow the existing safety data sheets.

Polymerisation times for Sebond pink:

Light curing device	Time
Spektra LED	1 min.
Spektra 2000 fast	3 min.
Spektra 2000	3 min.
DC-XS/Uni-XS	90 sec.
Hilite power	90 sec.
Spektramat	2 min.
GC Labolight LV-II/III	1 min.
Solidilite EX/V	1 min.

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If there is anything in this patient information leaflet that you do not understand, please contact our customer service department before using the product.

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