MINILUX

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

CONTENTS

1. GENERAL NOTES	10
2. DESCRIPTION OF THE UNIT AND ITS USE	11
3. TECHNICAL SPECIFICATIONS	11
4. IMPORTANT WARNINGS	11
5. INSTRUCTIONS FOR USE	12
6. MAINTENANCE	13
7. WARNINGS	13
8. WASTE DISPOSAL	13
9. MALFUNCTION CAUSES REMEDIES	14
10. SPARE PARTS	15

1. GENERAL NOTES

Please read this booklet carefully. It contains important information about safe installation, use and maintenance. Keep this booklet in a convenient place for further reference.

1.1 Check the general condition of the unit after removing the packaging. When in doubt do not operate the unit and contact a qualified service repair person authorised by **SIMED**.

1.2 Before connection, check to see that the identification plate indicates the correct voltage for the electrical supply system. Installation must be carried out in compliance with all applicable safety standards and in accordance with the manufacturer's instructions. The manufacturer will not be held responsible for damage to objects and/or harm to people caused by improper installation. The unit is safe only when connected correctly to ground via an appropriate device in accordance with safety standards. This fundamental requirement must be satisfied. If in doubt, request a thorough inspection of the system by a professionally qualified technician. The manufacturer will not be held responsible for damage caused by an improperly grounded system.

1.3 This unit should only be used as intended by the manufacturer. Any other use must be considered inappropriate and consequently dangerous. The manufacturer will not be held responsible for damage caused by improper, erroneous or unreasonable use.

1.4 Switch-off and unplug the unit from the electrical outlets before servicing.

1.5 Carefully follow cleaning instructions as in the manufacturer's manual.

1.6 In cases of breakdown or unit malfunction, turn the unit off and do not try to repair it. Contact only qualified and authorised SIMED agents. Any subsequent repairs made to the unit must be carried out by the manufacturer or by an authorised agent using original SIMED replacement parts only. Not following the above procedure could jeopardise the safety on the unit. Connection to the electrical distribution system must be made by using only plugs that conform to safety standards.

1.7 Check that the electrical supply system is adequate for the maximum power consumption/rating as indicated on the identification plate. When in doubt, contact a professionally qualified technician. The technician must check that the system cables are adequate for the current drawn by the equipment. The use of adapters, multiple plugs and extension cords is not advisable. If absolutely necessary, use devices that conform to safety standards, being careful not to exceed the limit values indicated on the device itself.


2. DESCRIPTION OF THE UNIT AND ITS USE

MINILUX is a unit for the electropolishing of non precious metal alloys. The power generator creates direct current which allows in the tank, filled with an appropriate solution (i.e. CROCOLIT SIMED), the electrolytic polishing of the element surface connected with the hook fastened to the central bar.

3. TECHNICAL SPECIFICATIONS


POWER SUPPLY

- Baked enamel coated steel housing.
- Voltage in: 230V 50Hz a.c. or 115V 60Hz a.c.
- Voltage out : 20 V d.c.
- Maximum power: 12A
- Overall dimensions (l x w x h) 22 x 20 x 21 cm
- Weight approx.: 6 kg

TANK

- Material: perspex
- Dimensions (l x w x h): 16 x 18,5 x 17 cm.
- Volume approx. 2,5 litres
- Voltage at the contacts can be regulated up to 24 V d.c.
- Copper cathodes for polishing CoCr-alloys or tungsten cathodes for polishing gold (optional)
- Titanium hook for hanging CoCr-frameworks or tungsten hook for hanging gold parts (optional)
- Peso ca.: 1,5 kg

4. IMPORTANT WARNINGS

4.1 Avoid any contact of the electropolishing liquid with water or humidity. Therefore avoid to immerse wet works into the bath or to add water. Furthermore it is good rule to work with the tank lid closed under a ventilating hood.

4.2 Once the electropolishing procedure is finished, bring the main switch (5) in "0"-position before removing the object from the tank.

4.3 Make sure that the hook, the anode bar and the object are in good contact in order to avoid sparking with overheating and local oxidation that might impair the electric current.

4.4 Check that all electric contacts are clean, free from oxidation and operationally effective.

4.5 In case of a short circuit or power surge, the automatic safety system will cut off the current immediately. In order to start up the unit again, press the reset button (20).


5. INSTRUCTIONS FOR USE

Connections

5.1 Check that the voltage of the main supply system is the same as required by the unit. Remove the shield (17)

5.2 Connect the tank cables (2) to the power supply terminals of the generator (1). Black cable (15) to the negative black terminal (18) and red cable (16) to the positive red terminal (21). Replace the shield (17).

5.3 Insert the power plug (19) into the power outlet.

Operating instructions

5.4 Fill the tank (2) 3/4 full with the electropolishing solution. The recommended electropolishing solution by Simed for CoCr-alloys is CROCOLIT.

5.5 Fix the framework on the hook (9) and hang it on the centre bar (12). Take care not to let the object touch the cathodes, and centre the object in the tank.

Attention: Clean perfectly the object and don't touch it with the naked hands before immersing it.

5.6 Turn on the power supply by putting the switch (5) to the "I" position. Set the timer (7) for the necessary time, about 8-10 minutes depending on the alloy, temperature and current intensity. More time is needed for a cold solution and less for a warm solution.

5.7 When the solution is cold, set the voltage regulator (4) to maximum (while the solution is cold only a little part of the current can pass through it) and little by little, as the solution heats up, reduce the power to keep the amperage between 3 and 5 A.


Attention: In case the object to be polished is small it could happen that the amperage doesn't raise (s.a. 9.2.3). In this case check if there is a formation of small bubbles near the cathodes to confirm that current flows.

5.8 The ideal temperature for good polishing results varies between 25 and 45/50 °C. Higher temperatures can cause the object to streak and stain and the hook to corrode. In such cases, let the solution cool down or dip the object into a tank with cold solution.

5.9 To speed up polishing when the solution is below 25 °C, dip a detective piece of an old framework into the solution by hanging it on the hook. Set the timer on approx. 15 minutes and set the voltage regulator to maximum. In this way you can preheat the solution.

5.10 If two frameworks are to be polished simultaneously, position them so that their surfaces are parallel to the cathodes and they do not touch or cover each other.

5.11 To improve the polishing properties and the penetrating power of the liquid as well as to reduce streaking effects and to minimise local overheating while maintaining an even temperature in the tank an electromagnetic stirrer can be used by placing it under the tank thereby acting as a rotating stirrer (optional).

6. MAINTENANCE

6.1 Minilux doesn't need any particular maintenance. Clean periodically the power supply unit, after disconnecting it from the main supply, with a dry cloth. If necessary use a not aggressive grease remover. When the electrolytic liquid must be changed clean the tank making attention to completely dry it before adding the new electrolytic liquid. Even few drops of water could damage the electrolytic solution.

7. WARNINGS

7.1 Absolutely avoid any contact between the central anode bar and the lateral cathodes. If the current is turned on while they are connected, it would cause a short-circuit with extreme sparking and severe overheating at the points of contact. Damage or breakdown of various power supply components would be the consequence.

7.2 Avoid touching the cathode plate with the hook or the framework.

7.3 It is good practice to first place the object to be polished in the tank and then turn the power on.

7.4 Do not lean over the tank while the equipment is in operation to avoid inhaling harmful fumes and being splashed with the solution. It is a good rule to work under a ventilating hood with the tank lid (8) closed.

7.5 If the solution spills, rinse immediately with plenty of water and then soak up the liquid using rags and sawdust.

7.6 The use of protective glasses and rubber gloves is recommended.

7.7 If the solution comes in contact with eyes or skin, rinse immediately with plenty of water. It is also possible to use a 20% sodium bicarbonate solution. If discomfort occurs seek medical advice.

8. WASTE DISPOSAL

8.1 The used solution is a harmful waste. Do NOT allow the product to be released into the environment. Follow the local regulations for the proper method of disposal.

This symbol means that the equipment, at the end of its life, must be disposed separately from the usual waste, following the local laws.

9. MALFUNCTION CAUSES-REMEDIES

9.1 Unit does not start up.

- 9.1.1 Check that the plug is inserted correctly and the cable is not damaged. Have the outlet checked by a qualified person to be sure that the voltage measures 220 Volt.
- 9.1.2 Check that the main switch (5) is functioning.
- 9.1.3 Check that the fuse is not burned out.

9.2 Too little or no electric current passes through the tank



- 9.2.1 Check if there are bubbles near the cathodes. If yes, current is passing through otherwise there is no current at all.
- 9.2.2 If current doesn't pass check if all the contacts are ok, that the main switch (5) and the light (6) are lighting, the timer (7) is started and the thermo-magnetic switch (20) is set. If everything is ok and still no current is passing contact SIMED or a qualified service repair person authorised by SIMED.
- 9.2.3 If current passes but the indication of the Amperemeter (3) doesn't rise, this could depend from the metal piece to be polished. If the piece to be polished is of reduced dimension the amperes passing will rise very slowly and could even not reach 3 Ampere, but the electropolishing effect is still working. In this case it is necessary to estimate at sight the polishing level or it can be tested with an external tester that the current passing is about 13-14 Volts. Otherwise check if all the contacts are ok, that the solution is getting warm (temperature 25°-30° C), that the main switch (5) and the light (6) are lighting and that the timer (7) and the power regulator (4) are working. e che il regolatore (4) ed il timer (7) funzionino. If necessary contact SIMED or a qualified service repair person authorised by SIMED.

9.3 The framework is not adequately polished.

- 9.3.1 Electrolytic solution too old or damaged (f.e. by water or humidity).
- 9.3.2 Solution too cold or too hot.
- 9.3.3 Solution not suitable for the metal alloy to be polished.
- 9.3.4 Surface of the framework not adequately prepared; presence of oxydation, weldings, insufficient sandblasting, fingerprints
- 9.3.5 Presence of excessive cavities where the current cannot reach the bottom. In case of big cavities (f.e. very arched frameworks) it is possible to try to solve the inconvenient by adding a thick copper thread to the cathode in order to approach the more distant point of the arched surface.
ATTENTION: Do absolutely NOT touch the hook, neither the framework nor other metal elements in contact with the hook (s.a. 7.1.)
- 9.3.6 To improve the polishing properties and the penetrating power of the liquid as well as to reduce streaking effects and to minimise local overheating while maintaining an even temperature in the tank an electromagnetic stirrer can be used by placing it under the tank thereby acting as a rotating stirrer.

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