

LASCOD

laboratory use



Progressing research and innovation



Plaster model has to reproduce with the highest fidelity every detail of the oral cavity detected by impression material. Progressing research in new products, accurate selection of raw materials and innovating working techniques are our daily business, in order to provide the perfect mix of quality, performance and time management. This allowed us to create the new Singletypo 4 Light Grey CAD SYSTEM plaster type IV extra-hard which widens the range of our plasters, well known by their high quality standards.

thixotropy and flowability

LASCOD plasters are characterized by excellent thixotropic properties. The product consistency is always at the desired level after mixing, allowing the perfect reproduction of details.

Excellent product flowability on the impression surface guarantees an accurate transfer of details.

Optimal results are achieved by pouring the material slowly in small volumes. The above mentioned approach favors the elimination of air-bubbles.




fast setting time and accuracy



Micronized and fine powder allows the accurate and detailed transfer of the oral cavity data recorded on the impression material to your model stone. An accurate model is your secret for a very good prosthesis.

Our plasters/stone are characterized by great strength and surface hardness, strong edges even in thin layers and no splinters during cutting or refining. No breakage risks when separating models from impressions.

CAD-CAM compatibility



Physical properties, specific color and innovative formulation make **Singletypo 4 Light Grey**, **Singletypo 4 Golden Brown**, **Kromotipo 4**, the perfect stone for CAD-CAM reading optical, laser and tactile.

enduring quality and long shelf life



Accurate selection of raw materials is our secret for supplying you with consistent quality for every batch of material.

Hermetic sealing of our buckets grant a long shelf life for the product without altering its properties. Lids of different colours make different kind of plasters easily noticeable.

managing and saving time

KTP090



With our working time regulator TIME OUT you can mix a larger quantity of plaster and cast several models at the same time. Extending it to your liking, LASCOD plasters/stones working time will not alter our products physical properties and performance in any way.

other Lascod laboratory products



EGX370

EGX070



EGX390

EGX090



XLG070



EGS160

EGS001



ALC100



The whole range of laboratory instruments can be found on Zeffiro General Catalogue or onto the website www.lascod.it at 'instruments' page.



CAD SYSTEM

TYPE IV EXTRAHARD

Master models, crown and bridges stumps, inlay/on-lay, metal frameworks, implants, antagonist.

LIGHT GREY

Mixing ratio (powder/water)	100/22
Imbibition time	20"
Manual mixing time	60"
Mechanical mixing time	30"
Pouring time	7'
Setting time	14'
Extraction time	30'
Expansion setting	0,09%

Resistance to compression	103 MPa	
	14.927 Psi	
1.051 Kg/cm ²		
Rockwell hardness	95HRI	
Packaging	25 Kg	TXG425
	6 Kg	TXG406
	10x1,250 Kg	-



CAD SYSTEM

TYPE IV EXTRAHARD

Master models, crown and bridges stumps, inlay/on-lay, metal frameworks, implants, antagonist.

GOLDEN BROWN

Mixing ratio (powder/water)	100/20
Imbibition time	20"
Manual mixing time	60"
Mechanical mixing time	30"
Pouring time	3'
Setting time	7'
Extraction time	30'
Expansion setting	0,08%

Resistance to compression	84 MPa	
	12.174 Psi	
860 Kg/cm ²		
Rockwell hardness	90HRI	
Packaging	25 Kg	TXE425
	6 Kg	TXE406
	10x1,250 Kg	TXE412




CAD SYSTEM

TYPE IV EXTRAHARD

Master models, crown and bridges stumps, inlay/on-lay, metal frameworks, implants, antagonist.

CHROMATIC - PASTEL



Mixing ratio (powder/water)	100/22
Imbibition time	20"
Manual mixing time	60"
Mechanical mixing time	30"
Pouring time	3'
Setting time	7'
Extraction time	30'
Expansion setting	0,08%

Resistance to compression	81 MPa	
	11.740 Psi	
828 Kg/cm ²		
Rockwell hardness	85HRI	
Packaging	25 Kg	KTP096
	6 Kg	KTP098
	10x1,250 Kg	KTP412



TYPE IV EXTRAHARD

Master models with complete and partial prostheses, metal frameworks, implants, antagonist.

PINK

Mixing ratio (powder/water)	100/25
Imbibition time	20"
Manual mixing time	60"
Mechanical mixing time	30"
Pouring time	3'
Setting time	7'
Extraction time	30'
Expansion setting	0,08%

Resistance to compression	70 MPa	
	10.144 Psi	
713 Kg/cm ²		
Rockwell hardness	75HRI	
Packaging	25 Kg	TPM425
	6 Kg	-
	10x1,250 Kg	TPM412



TYPE III HARD

Antagonist models, Master models with complete and partial prostheses, metal frameworks, study and preliminary models.

GREEN

YELLOW

Mixing ratio (powder/water)	100/29
Imbibition time	20"
Manual mixing time	60"
Mechanical mixing time	30"
Pouring time	3'
Setting time	7'
Extraction time	30'
Expansion setting	0,09%


Resistance to compression	66 MPa	
	9.565 Psi	
670 Kg/cm ²		
Rockwell hardness	70 HRI	
Packaging (YELLOW)	25 Kg	TGI325
	6 Kg	TGI306
	10x1,250 Kg	TGI312
Packaging (GREEN)	25 Kg	TGV325
	6 Kg	-
	10x1,250 Kg	TGV312



TYPE III HARD

Antagonist models, Master models with complete and partial prostheses, metal frameworks, study and preliminary models.

CHROMATIC - LIGHT BLUE



Mixing ratio (powder/water)	100/29
Imbibition time	20"
Manual mixing time	60"
Mechanical mixing time	30"
Pouring time	3'
Setting time	7'
Extraction time	30'
Expansion setting	0,09%

Resistance to compression	66 MPa
	9.565 Psi
Rockwell hardness	70HRI
Packaging	25 Kg KTP095
	6 Kg KTP097
	10x1,250 Kg KTP312



TYPE IV EXTRAHARD

Orthodontic, demonstration and antagonist.

BRILLIANT WHITE

Mixing ratio (powder/water)	100/22
Imbibition time	20"
Manual mixing time	60"
Mechanical mixing time	30"
Pouring time	3'
Setting time	7'
Extraction time	30'
Expansion setting	0,08%

Resistance to compression	103 MPa
	14.927 Psi
Rockwell hardness	95HRI
Packaging	25 Kg TOR425
	6 Kg -
	10x1,250 Kg TOR412



TYPE III HARD

Orthodontic, demonstration and antagonist.

BRILLIANT WHITE

Mixing ratio (powder/water)	100/25
Imbibition time	20"
Manual mixing time	60"
Mechanical mixing time	30"
Pouring time	3'
Setting time	7'
Extraction time	30'
Expansion setting	0,09%

Resistance to compression	60 MPa
	8.700 Psi
Rockwell hardness	60HRI
Packaging	25 Kg TOR325
	6 Kg TOR306
	10x1,250 Kg TOR312



TYPE III HARD

Setting on articulator, metal framework, fast setting, low expansion, high adhesiveness.

WHITE

Mixing ratio (powder/water)	100/30
Imbibition time	20"
Manual mixing time	60"
Mechanical mixing time	30"
Pouring time	3'
Setting time	4'
Extraction time	30'
Expansion setting	0,05%

Resistance to compression	48 MPa
	6.956 Psi
Rockwell hardness	-
Packaging	20 Kg TRT220
	6 Kg -
	10x1,250 Kg TRT212



TYPE II SOFT

Flasking and creation of full or partially removable prosthesis, setting on articulator, metal frameworks, study models.

WHITE

Mixing ratio (powder/water)	100/43
Imbibition time	20"
Manual mixing time	60"
Mechanical mixing time	30"
Pouring time	3'
Setting time	7'
Extraction time	30'
Expansion setting	0,09%

Resistance to compression	29 MPa
	4.203 Psi
Rockwell hardness	-
Packaging	20 Kg TMF220
	6 Kg -
	10x1,250 Kg -

trouble shooting guide

F.A.Q.

Why isn't my plaster/stone setting quickly enough?

- Make sure to thoroughly clean the impression from residual traces of blood and saliva.
- Make sure that there is no left-over water deposit on the impression after rinsing.
- Your impression material may not be compatible with your plaster.
- Make sure to use the water/powder ratio suggested by the manufacturer. Did you use too much water?
- Make sure to stir the plaster powder in the original packaging or in the drawer before use.
- Make sure to follow the manufacturer's instruction for manual and/or mechanical mixing.
- Make sure to store your plaster container away from direct sunlight, heating sources and in a cool and dry environment.
- Avoid the use of soap or liquid detergent to clean the mixing bowl.
- Try using the mixing water at room temperature.
- Plaster/stone powder should not absorb humidity before use. Make sure to close the container lid tightly.

Why is my plaster/stone setting too quickly?

- Make sure that mixing bowl and spatulas used for manual or mechanical mixing are perfectly clean.
- Make sure to use the water/powder ratio suggested by the manufacturer. Did you use less water?
- Make sure to follow the manufacturer's instruction for manual and or mechanical mixing.
- Make sure to use water at room temperature and not to work at extreme temperatures.
- Make sure that your tap water is not too hard. Did you use the model trimmer waste water?
- Avoid extending the vibrating time while pouring your plaster/stone on the impression.
- Have you used excessive amounts of salt or setting accelerators?
- Maybe you have submerged the impression in potassium sulfate base solution for too long.
- Make sure to store your plaster container away from direct sunlight, heating sources and in a cool and dry environment.
- Your plaster/stone contains hardened residuals.

Why isn't my model surface homogenous and accurate?

- Make sure that mixing bowl and spatulas used for manual or mechanical mixing are perfectly clean.
- Make sure to store your plaster container with its lid tightly closed.
- Make sure to eliminate air-bubbles during manual mixing and/

or check the your mixer vacuum is working correctly.

- When pouring the powder in the mixing bowl, try to do it slowly and in small quantities to facilitate the elimination of air-bubbles.
- Make sure to wait at least 30 minutes before separating model and impression.
- Make sure to thoroughly clean the impression from residual traces of blood and saliva.
- Avoid using excessive power on the vibrating device.
- Make sure to follow the manufacturer's instruction for manual and or mechanical mixing. If your model show stripes of different shades, you may want to set your vacuum mechanical mixing device on a longer mixing time.
- After trimming, make sure to eliminate left over waste on the model with a soft brush.
- Make sure not to use excessive amounts of salt or setting accelerators as well as model trimmer waste water. Sodium chloride will increase expansion and calcium sulphate will stain your model surface.
- When your plaster loses its shine on surface, the working time is almost over. Avoid modelling for more than 1 minute.
- Your impression material may not be compatible with your plaster.
- Against rounded edges avoid an excessive use of the steam cleaner.
- Check that the impression received is sufficiently accurate.
- Try using mechanical mixing with vacuum.
- Try modeling the model without rushing.
- Avoid placing your model on paper while it is setting.
- Make sure that your model is not completely dry when you use the steam cleaner or before submerging it in boiling water.

Why isn't my plaster/stone hard enough?

- Avoid using excessive amounts of water.
- Avoid using excessive manual or mechanical mixing time.
- Avoid adding too much salt to water used in the mix.
- Avoid working on impression with residual traces of blood and saliva.
- Avoid placing the mixing bowl on the vibrator.
- Wait at least 30 minutes before separating model and impression.
- Stir the plaster powder in your drawer before use.
- Use compatible products.

The photographic images here reproduced are purely indicative and are not necessarily identical to the actual products.



For further
information
activate the
QR Code reader
on your
smartphone.

