discovery
pearl

Gebrauchsanweisung
Instructions for use
Mode d'emploi
Modo de empleo
Modalità d'uso

D
DENTAURUM

- Informationen zu Produkten finden Sie unter www.dentaurum.com
- For more information on our products, please visit www.dentaurum.com
- Vous trouverez toutes les informations sur nos produits sur www.dentaurum.com
- Descubra nuestros productos en www.dentaurum.com
- Informazioni su prodotti sono disponibili nel sito www.dentaurum.com

Stand der Information | Date of information | Mise à jour |
Fecha de la información | Data dell'informazione: 04/18

989-301-04 Printed by Dentaurum Germany 0418/BRT

D
DENTAURUM

Turnstr. 31 | 75228 Ispringen | Germany | Tel. +49 72 31 / 803-0 | Fax +49 72 31 / 803-295
www.dentaurum.com | info@dentaurum.com

DE

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EN

Dear customer,

Thank you for choosing a quality product from Dentaurum.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit.

Instructions for use cannot describe every eventuality and possible application. Should you have any questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the Internet at www.dentaurum.com even though you may frequently use the same product.

1. Manufacturer

Dentaurum GmbH & Co. KG | Turnstr. 31 | 75228 Ispringen | Germany

2. General description

discovery® pearl is a ceramic bracket manufactured using the CIM process (Ceramic Injection Molding). The bracket is made from highly pure polycrystalline aluminum oxide.

3. Intended use

For the orthodontic treatment of malocclusions.

4. Processing

4.1 Informing the patient

Inform the patient about the possible consequences and risks of bracket treatment, in particular during debonding. This applies especially in the case of teeth with enamel that already exhibits signs of initial damage (cracks, crazing, splitting, etc.).

4.2 Conditioning teeth

Follow the instructions from the bracket adhesive manufacturer for conditioning the teeth.

4.3 Selection of bracket adhesives

discovery® pearl brackets can be bonded using all commercially available bracket adhesives. We recommend CONTEC lc (REF 163-111-00), ConTec Go! (REF 163-321-00) and SmartBond® (REF 163-506-00).

4.4 Handling brackets

Grip the discovery® pearl brackets mesially-distally using clamping tweezers (e.g. REF 025-277-00) to remove them from the packaging and handle according to the following steps. The base of the brackets should not be touched with your fingers or contaminated in any way, as this could otherwise detrimentally affect the bond.

The brackets can be identified and matched to the correct tooth by means of the FDI marking (e.g. 21) on their base.

4.5 Applying bracket adhesive

Prepare the teeth and the bracket adhesive according to the manufacturer's instructions. Apply a small amount of adhesive to the center of the bracket base, enough to wet the entire bracket base. This guarantees an optimum bond.

4.6 Indirect bonding technique

discovery® pearl brackets are bonded on the plaster or the plastic model with the bracket adhesive. The adhesive layer may not be removed afterwards! Indirect bonding techniques using a sugar solution or other water-soluble lab adhesives that later have to be completely removed are not suitable. This would destroy the mechanical retention of discovery® pearl and a reliable bracket fit could no longer be guaranteed.

Thermoformed splints or silicone transfers, e.g. LinguSil (REF 164-000-00) can be used as a transfer matrix.

4.7 Placing and aligning brackets

After applying the adhesive, place, align and press the bracket in the usual way onto the tooth until the adhesive exudes from the sides.

Avoid contact with the antagonist due to the risk of abrasion and/or premature bracket loss. If necessary, a palatal bite raise can be used e.g. Bite-Bumper™ (REF 705-000-00).

Remove completely any excess adhesive before it has cured. This is particularly important to ensure easy debonding at the end of treatment.

4.8 Ligating the archwire

We recommend using Tensic®, rematitan® LITE or the translucent pearl archwires at the beginning of treatment. **Do not use force to ligate the archwires. Do not torque the archwires when they are ligated.**

Once the adhesive has cured (curing time according to the adhesive manufacturer's instructions), ligating can be completed in the usual way using wire or plastic ligatures (e.g. Dentalastics® pearl REF 774-563-00).

Note: Careful ligation avoids damage to the brackets!

4.9 Bracket debonding

First, remove the archwire. For debonding the discovery® pearl brackets we recommend the special discovery® pearl debonding instrument (REF 019-001-00) which grips the body of the bracket directly at the base. The bracket is gently removed from the bracket adhesive by slightly tipping it sideways in the direction of the slot and archwire.

Alternatively, the straight bracket removal pliers (REF 004-346-00) can be used for debonding the discovery® pearl brackets. Place the blades of the bracket removal pliers mesio-distally on the bracket base, **directly on the surface of the tooth**. The bracket is gently removed from the tooth by slightly twisting the pliers.

Note:

Despite their small, delicate shape there is no risk of fracture with discovery® pearl brackets provided they are used correctly. Use of other instruments can result in bracket fracture. If the bracket fractures during incorrect debonding, we recommend carefully removing the fractured segments and the adhesive from the enamel using bracket debonding pliers (e.g. REF 003-349-00) or a scaler (e.g. REF 027-349-00).

Use care when debonding the segments, as there is a risk of the enamel chipping. Use an aspirator and wear protective glasses.

5. Composition

The composition is included in the materials list, refer to the catalog or www.dentaurum.com.

6. Safety information

The product should not be used if there is a known allergic reaction to one or more of the material components.

7. Storage and shelf life

Store discovery® pearl brackets in closed packaging.

discovery® pearl brackets do not have a limited shelf life.

8. Information for single use products

discovery® pearl brackets are intended for single use only. Reconditioning of discovery® pearl brackets that have already been used (recycling) and their reuse on a patient is not permitted.

Micro cracks may result from loading during treatment and in particular bracket debonding that could cause fractures in the ceramic if reused!

9. Additional information

- Technical information such as torque, angulation, size etc. can be found on the relevant product label or in the catalog information.
- There are no particular precautions to be taken when disposing of discovery® pearl brackets after treatment.

10. Availability

For the complete range of discovery® pearl brackets, please see our catalog or www.dentaurum.com.

11. Quality information

Dentaurum ensures faultless quality of the products manufactured by us. These recommendations are based upon Dentaurum's own experiences. The user is responsible for the processing of the products. Responsibility for failures cannot be taken by Dentaurum as we have no influence on the processing on site.

12. Explanation of symbols used on the label

 Refer to the label. Additional information can be found at www.dentaurum.com (Explanation of symbols REF 989-313-00).

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