


Mixing tips

Apply the mixing tips on the cartridge and for locking turn the bayonet lock. It is recommended to dispose the initially mixed (pea-sized) amount of material.



Important Advice

When using double cartridges, the mixing tip is to remain onto the cartridge as sealing.

Form of delivery

			
Mixing ratio 1:1 Inner diameter 5.4 mm Colour transparent	Mixing ratio 1:1 Inner diameter 5.4 mm Colour pink	Mixing ratio 1:1 Inner diameter 6.5 mm Colour green	Mixing ratio 1:1 Inner diameter 3.2 mm Colour light blue
Product only for Dospur M4	Product Dospur evo/M5 StoneBite® StoneBite® scan StoneBite® scan plus Sifix	Product Fresh® clear Fresh® heavy body Fresh® mono Regofix MegaBase® Odontosil	Product AdiFilm LP GumQuick GumQuick scan GumQuick implant
40 pieces D-4960 100 pieces D-4960-1	40 pieces D-4970-1 100 pieces D-49701-1	40 pieces D-4970-2 100 pieces D-49701-2	40 pieces D-4970-3 100 pieces D-49701-3
			
Mixing ratio 4:1 Inner diameter 3.2 mm Colour light blue	Mixing ratio 1:1 Inner diameter 4.2 mm Colour yellow	Mixing ratio 1:1 Inner diameter 2.5 mm Colour brown	Mixing ratio 5:1 Inner diameter 6.0 mm Colour yellow
Product Fixtemp C&B	Product Fresh® light	Product Fixtemp cement	Product Fresh® 380 ml Zerisil® soft
40 pieces D-4970-4 100 pieces D-49701-4	40 pieces D-4970-5 100 pieces D-49701-5	40 pieces D-4970-6 100 pieces D-49701-6	40 pieces D-4981 100 pieces D-4982

