



AMANNGIRRBACH

DESIGNED TO
OUTPERFORM.

 **zolid gen-x**


THE UNIVERSAL ZIRCONIA THAT CHANGES EVERYTHING.



 **AMANNGIRRBACH AG**
Herrschaftswiesen 11 6842 Koblach | Austria

UNIVERSAL EXCELLENCE

The new Zolid Gen-X is a multilayer, highly translucent zirconia and is a true all-rounder in terms of "universal use". Due to its high strength and natural translucency, this material can be used for virtually every indication. Even the massive structures of implant-supported restorations with a gingival component appear completely natural. Zolid Gen-X is available in the 16 classic VITA tooth shades as well as two Bleach shades. This makes the fabrication of restorations easy and efficient for every laboratory.






THE HIGHLIGHTS AT A GLANCE

- Best choice:** reduces the complexity of choosing the right material.
- Naturally beautiful:** the flowing color and translucency gradient is a perfect imitation of nature.
- Maximum flexibility:** due to a flexural strength of over 1,000 MPa, virtually all types of indications can be realized.
- New species:** Zolid Gen-X combines new properties with all the proven benefits of the Zolid HT+ family in a single material.

OPTIMIZED WORKFLOW

Zolid Gen-X can be integrated easily into the daily laboratory routines and is perfectly adapted to the workflow in daily laboratory work. This relates to milling processes, standard sintering cycles in the Zolid DNA workflow, as well as easy finalization with the 16 A-D VITA shades and Ceramill Stain & Glaze. This greatly facilitates the fabrication process for the laboratory.

<p>Perfect handling</p>  <p>Perfectly matched to the dental milling process without making any compromises.</p>	<p>Standard sintering cycle</p>  <p>Sintering in the Zolid DNA workflow at 1450 °C/2 h</p>	<p>Easy fabrication</p>  <p>Easy finalization with 16 A-D VITA shades and Ceramill Stain & Glaze</p>
--	---	---



"Zolid Gen-X is a zirconia which combines both the worlds of strength and esthetics. It allows me to plan and execute small to complex reconstructions with full control over the restoration."

CDT Alexander Wuensche
Zahntechnik Inc. | Miami FL, USA




THE PEAK OF QUALITY

Zolid DNA blanks are manufactured 100 percent in-house at the Amann Girschbach manufacturing site Dentustry One in Austria. Only the best raw materials are processed to blanks using the strictest test methods and specified, certified processes. This ensures that all material parameters are perfectly matched to the needs of dentistry. Fully embedded into the Ceramill CAD/CAM workflow, this results in an efficient and safe system with tailor-made solutions.




<p>100%</p>  <p>MANUFACTURED IN-HOUSE</p>	<p>10 MIO UNITS</p>  <p>CLINICALLY PROVEN</p>	<p></p> <p>OPEN FOR ALL SYSTEMS</p>	<p></p> <p>TOSOH</p>
---	---	--	---

 100% Made in Austria

ZOLID DNA PORTFOLIO ZIRCONIA FOR ALL REQUIREMENTS

zolid gen-x



NEW
16 A-D VITA SHADES

"One for all" material with best mechanical properties and natural appearance, suitable for a wide range of indications


zolid fx multilayer



16 A-D VITA SHADES

Unbeatable zirconia for esthetic demands and a natural appearance in the anterior region


zolid ht+ preshades



16 A-D VITA SHADES

The best option for highly efficient multi-unit restorations

zi white



For fully veneered frames with individualized esthetics

zolid ht+ white



Ideal for REAX restorations with an individual touch

zolid fx white



Best solution for individually colored, highly esthetic anterior restorations

TRANSLUCENCY →

EFFICIENCY →

ORDERING AND PRODUCT INFORMATION

Further product and ordering information on Zolid Gen-X and the Zolid DNA Generation can be found in the following media:

- Digital Lab Equipment (Art. No. 9FMD19),
- Ceramill Material - Materials for CAD/CAM processing (Art. No. 976023),
- Zolid DNA Generation processing techniques (only electronically).

All are for as download on our website.



www.amanngirrbach.com/zolid-gen-x/